

PRODUCT
CATALOGUE

bpositive

MAKING
A **POSITIVE**
IMPACT
TOGETHER



We can all be heroes by helping others with
buy1help1





IMAGINE A WORLD WHERE EVERYONE IS HEALTHY AND POSITIVE.

Imagine a company that wants to work with you to make this happen. To make a positive impact together.

You are making a positive impact each time you use our products. Firstly, by healing your patients. Secondly, every product purchased by you helps someone else through our buy1help1 initiative.



To learn more, and choose which charity to support, visit our website or scan the QR Code.

SUMMARY OF PRODUCT RANGE



Absorbent Dressings (Standard & Lite)



Alcohol Wipes



Bandaging



Combine Dressings / Gamgee



Plasters (Spot & Strip)



Gauze Swabs - Non Woven



Gauze Swabs - Woven



Gloves



Masks



Soft & Flexible Dressings & Roll



Water & Bacteria Proof Dressings & Roll



Emesis Bag



Hydrocolloid Dressing






Closure Strips



Woundcare Retail Packs

Boo Boo Plasters Minions Edition



Flexible



Protect



Comfortable

© Universal City Studios LLC. All Rights Reserved.



BOO BOO PLASTERS			
Code	Description	PCS/BX	BX/CRT
011320WRA	Boo Boo Plaster Strips for Adults (BX20)	20	48
011340KRA	Boo Boo Plaster Strips for Kids (BX40)	40	48
011320KRA	Boo Boo Plaster Strips for Kids (BX20)	20	48
011310KRA	Boo Boo Plaster Strips for Kids (BX10)	10	240
010604RA	Boo Boo Special Dressing 4 x 5cm (PK5 retail)	5	180



Painless Removal



Helps Protect Against Germs



Absorbs And Cushions



Hospital Grade

ABSORBENT NON-ADHERENT Dressing

For cuts, grazes and minor burns

This ABSORBENT NON-ADHERENT DRESSING helps wounds to heal with a specially designed absorbent pad which has a soft cushioning effect to protect the wound. The pad absorbs excess wound fluid which is important for wound healing. We have added an extra layer to make it easier and less painful when changing the dressing. Keep the dressing in place with a bandage, adhesive tape or waterproof film dressing.

BPOSITIVE N/ADH ABSORBENT DRS

Code	Description / sizes	PCS/BX	BX/CRT
010105BA	N/ADH ABSORBENT DRS 5 X 5CM	100	20
010110BA	N/ADH ABSORBENT DRS 10 x 10CM	100	8
010120BA	N/ADH ABSORBENT DRS 10 X 20CM	50	6
010105RA	N/ADH ABSORBENT DRS 5 X 5CM (Retail)	4	270
010110RA	N/ADH ABSORBENT DRS 10 X 10CM (Retail)	3	120



Painless Removal



Helps Protect Against Germs



Absorbs And Cushions



Hospital Grade

ABSORBENT NON-ADHERENT Lite Dressing

For cuts, grazes and minor burns

This LITE ABSORBENT NON-ADHERENT DRESSING is the same construct as our standard version but has approximately 10-15% less cotton. A good step-down where the wound has less discharge.

BPOSITIVE LITE N/ADH ABSORBENT DRS

Code	Description / sizes	PCS/BX	BX/CRT
010205BA	LITE N/ADH ABSORBENT DRS 7.5 X 5CM	100	12
010210BA	LITE N/ADH ABSORBENT DRS 7.5 X 10CM	50	16
010220BA	LITE N/ADH ABSORBENT DRS 7.5 X 20CM	50	12

ALCOHOL SWABS

Keeps you clean

The primary use of our alcohol swab is to clean the skin before injections, blood tests, or minor surgical procedures. The alcohol in the pad helps to kill or reduce the number of microorganisms. The pad is individually packaged in a aluminium foil wrapper to maintain its cleanliness and effectiveness.



BPOSITIVE LITE N/ALCOHOL WIPES

Code	Description / sizes	PCS/BX	BX/CRT
083670BA	70% ALCOHOL SWABS (Bx200)	200	50
083670RA	70% ALCOHOL SWABS (Bx100)	100	100






-  Gentle to Skin
-  Soft and Conformable
-  High Stretch
-  Hospital Grade



-  Gentle to Skin
-  Soft and Conformable
-  High Stretch
-  Hospital Grade



-  Moderate Compression
-  Gentle to Skin
-  Hospital Grade

BANDAGE Conforming Gauze

For securing dressings and light compression

This bandage is highly elastic, soft and conformable to allow joint movement. The lightweight construction provides good breathability and air circulation. Good for retaining dressings, light support and compression. Non-sterile and latex free.

-2m unstretched / ~3.8m stretched.-

BPOSITIVE CONF GZ BDG			
Code	Description / sizes	PCS/BX	BX/CRT
011103BA	CONF GZ BDG 2.5CM X 2M	12	160
011105BA	CONF GZ BDG 5CM X 2M	12	84
011175BA	CONF GZ BDG 7.5XM X 2M	12	60
011110BA	CONF GZ BDG 10CM X 2M	12	48
011115BA	CONF GZ BDG 15CM X 2M	12	36

BANDAGE Elastic Crepe

For securing dressings and light compression

This bandage is soft and conformable to allow joint movement. Good for retaining dressings, light support and compression. Non-sterile and latex free.

-1.8m unstretched / ~3.8m stretched.-

BPOSITIVE CREPE ELASTIC BDG			
Code	Description / sizes	PCS/BX	BX/CRT
0111003BA	CREPE ELASTIC BDG 2.5CM X 1.8M	12	81
0111005BA	CREPE ELASTIC BDG 5CM X 1.8M	12	54
0111075BA	CREPE ELASTIC BDG 7.5CM X 1.8M	12	36
0111010BA	CREPE ELASTIC BDG 10CM X 1.8M	12	27
0111015BA	CREPE ELASTIC BDG 15CM X 1.8M	12	18

BANDAGE Heavy Crepe

For securing dressings and supporting strains

This bandage is ideal for securing a wound dressing, supporting strains, and providing compression to alleviate swelling. Non-sterile and latex free.

-1.8m unstretched / ~3.4m stretched.-

BPOSITIVE CREPE HEAVY BDG			
Code	Description / sizes	PCS/BX	BX/CRT
011205BA	CREPE HVY BDG 5CM X 1.8M	12	42
011275BA	CREPE HVY BDG 7.5CM X 1.8M	12	30
011210BA	CREPE HVY BDG 10CM X 1.8M	12	24
011215BA	CREPE HVY BDG 15CM X 1.8M	12	12



COMBINE DRESSINGS / GAMGEE

This dressing consisting of a highly absorbent thick layer of fleece enclosed in a soft and conformable non-woven fabric, folded over to form a flat tube. Comes in Sterile & Non-Sterile versions.

This soft and absorbent dressing is used for highly exudating and/or bleeding wounds as a primary or secondary dressing.

BPOSITIVE CMBN DRE			
Code	Description / sizes	PCS/BX	BX/CRT
013410BA	Sterile Non-Woven Combine 10x12cm	50	10
013420BA	Sterile Non-Woven Combine 20x22cm	50	6
013422BA	Sterile Non-Woven Combine 10x22cm	50	3
013510BA	Non-Woven Combine Roll 10cmx10m	1	24
013520BA	Non-Woven Combine Roll 20cmx10m	1	12



DRESSING Spots & Strips

For protecting and covering small wounds, cuts, punctures or skin conditions such as warts.

BPOSITIVE DRESSING STRIPS / SPOTS			
Code	Description	PCS/BX	BX/CRT
0113100BA	PLASTER STRIPS (100)	100	24
011340RA	PLASTER STRIPS (Bx40)	40	48
011320RA	PLASTER STRIPS (Bx20)	20	48
0119100BA	FABRIC STRIPS (100)	100	24
011940RA	FABRIC STRIPS (Bx40)	40	48
011920RA	FABRIC STRIPS (Bx20)	20	48
0114100BA	PLASTER SPOTS (100)	100	24
01141000BA	PLASTER SPOTS (1000)	1000	6
011310RA	Plaster Strips (Bx10)	10	240
011320WRA	Special Edition Plaster Strips (Bx20)	20	120
011320ARA	Boo Boo Plaster Strips for Adults (BX20)	20	48
010320ARA	Boo Boo Transparent Plaster Strips for Adults (BX20)	20	48
011340KRA	Boo Boo Plaster Strips for Kids (BX40)	40	48
011320KRA	Boo Boo Plaster Strips for Kids (BX20)	20	48
011310KRA	Boo Boo Plaster Strips for Kids (BX10)	10	240



Flexible



Protect



Comfortable



FEEL BETTER SOON - EMESIS BAG			
Code	Description	PCS/BX	BX/CRT
052800BA	EMESIS BAG	50	12



GAUZE SWABS - NON WOVEN

NON-STERILE GAUZE NW			
Code	Description / Sizes	PCS/BX	BX/CRT
011805NWA	NW GAUZE 5 x 5CM (4 PLY) 30gm	100	50
011875NWA	NW GAUZE 7.5 x 7.5CM (4 PLY) 30gm	100	50
011810NWA	NW GAUZE 10 x 10CM (4 PLY) 30gm	100	40
STERILE GAUZE NW			
011875NWSA	Sterile NW GAUZE 7.5 x 7.5CM (4 PLY) PK5 40gm	1	600
011875NWS3A	Sterile NW GAUZE 7.5 x 7.5CM (4 PLY) PK3 40gm	1	600



GAUZE SWABS - WOVEN

This soft and highly absorbent gauze can be used to clean a wound before it is dressed. Comes in Sterile & Non-sterile versions.

NON-STERILE GAUZE SQ			
Code	Description / Sizes	PCS/BX	BX/CRT
010805BA	GAUZE SQ 5CM 8PLY	100	200
010875BA	GAUZE SQ 7.5CM 8PLY	100	100
010810BA	GAUZE SQ 10CM 8PLY	100	100
STERILE GAUZE SQ			
010975DA	GAUZE SQ 7.5 X 7.5CM 5'S STRL D/PACK*	50	6
010975YA	GAUZE SQ 7.5 X 7.5CM 3'S STRL	1	500
010975XA	GAUZE SQ 7.5 X 7.5CM 5'S STRL	1	500
010910BA	GAUZE SQ 10 X 10CM 5'S STRL	1	500
010810BG8A	GAUZE SWABS 10CM (8 PLY) (GREEN)	100	60
010810BG16A	GAUZE SWABS 10CM (16 PLY) (GREEN)	100	40

*Dispenser Pack



BIODEGRADEABLE NITRILE GLOVES

Want to do your part to help the environment? Typical nitrile gloves can take over 300 years to biodegrade whereas ours take between 2-5 years. They are intended to be worn by healthcare personnel who conduct medical examinations, during diagnostic and therapeutic procedures, and also for laboratory work. The gloves can also be used in clinical, home, or industrial environments.

BIODEGRADEABLE NITRILE GLOVES			
Code	Description / Sizes	PCS/BX	BX/CRT
0732XSBA	BNP FREE GLOVES XS (3.0 mil) 3.5g		
0732SBA	BNP FREE GLOVES S (3.0 mil) 3.5g		
0732MBA	BNP FREE GLOVES M (3.0 mil) 3.5g		
0732LBA	BNP FREE GLOVES L (3.0 mil) 3.5g		
0732XLBA	BNP FREE GLOVES XL (3.0 mil) 3.5g		



- 
 Optimal Wound Healing
- 
 Protects wounds from infection
- 
 Waterproof for showering
- 
 Protects minor wounds



- 
 Optimal Wound Healing
- 
 Protects wounds from infection
- 
 Waterproof for showering
- 
 Protects minor wounds



SCIENTIFICALLY ADVANCED HEALING HYDROCOLLOID DRESSING

Our scientifically advanced dressing helps promote faster healing, prevent scarring and reduce pain. The top layer offers protection against water, dirt and germs. The intended purpose of our dressing is to provide an optimal condition for wound healing and create a barrier to protect wounds from bacteria and pathogens. It is designed for dry to low exuding wounds. This is not intended to be a primary dressing used on chronic wounds or wounds involving the dermis.

BPOSITIVE SAH HYDROCOLLOID DRS			
Code	Description / Sizes	PCS/BX	BX/CRT
061905RA	SAH HYDROCOLLOID DRESSING 5CM x 5CM (retail)	5	
061908RA	SAH HYDROCOLLOID DRESSING 7.5CM x 7.5CM (retail)	5	
061910RA	SAH HYDROCOLLOID DRESSING 10CM x 10CM (retail)	3	
061904RA	SAH HYDROCOLLOID OVAL DRESSING 4.4 x 6.9CM (retail)	5	
061905BZ	SAH HYDROCOLLOID DRESSING 5CM x 10CM	10	
061910BZ	SAH HYDROCOLLOID DRESSING 10CM x 10CM	10	
061905BY	SAH HYDROCOLLOID DRESSING 5CM x 20CM	10	

ADVANCED HEALING HYDROCOLLOID DRESSING (Plus)

Our advanced hydrocolloid dressing helps promote faster healing, prevent scarring and reduce pain. The intended purpose of our dressing is to provide an optimal condition for wound healing and create a barrier to protect wounds from bacteria and pathogens. It can be used in the management of low to moderate exuding wounds like venous leg ulcers, pressure ulcers, superficial partial – thickness burns, postoperative wounds, skin abrasions and small wounds. It can also be used to rehydrate necrotic tissue, which is then removed by autolysis.

BPOSITIVE AH HYDROCOLLOID DRS (Plus)			
Code	Description / Sizes	PCS/BX	BX/CRT
061905BA	AH HYDROCOLLOID DRESSING (Plus) 5CM x 5CM	10	
061910BA	AH HYDROCOLLOID DRESSING (Plus) 10CM x 10CM	10	

MASKS

The Level 2 and 3 Surgical Face Mask is designed for single-use in healthcare environments where individuals may be exposed to the risk of blood or fluid splashes. These masks effectively block the passage of blood or bodily fluids through their multiple layers. Furthermore, they offer both comfort and breathability, featuring two ear loops that securely hold the mask in position.

NON-STERILE SURGICAL MASKS			
Code	Description / Sizes	PCS/BX	BX/CRT
011550BA	L1 Ear Loop 3 Layers	50	10
011550XA	L2 Ear Loop 3 Layers	50	10
011550YA	L3 Ear Loop 3 Layers	50	10
031550ZA	L2 Ear Loop 4 Layers (Australian Made)	50	27



- 
 For lacerations & surgical incisions
- 
 Support sutures and staples
- 
 Close your wound faster
- 
 Reduce scarring

SECURE WOUND CLOSURE STRIPS

The intended use of this product is to help lacerations and surgical incisions heal by bringing the edges of the wound together to close. Our product may be used in conjunction with skin sutures and staples.

BPOSITIVE SW CLOSURE STRIPS			
Code	Description / Sizes	PCS/BX	BX/CRT
022206BA	SW CLOSURE STRIPS 6MM x 38MM	50	
022206BZ	SW CLOSURE STRIPS 6MM x 75MM	50	
022212BA	SW CLOSURE STRIPS 12MM x 100MM	50	
022206RA	SW CLOSURE STRIPS 6MM x 38MM (Retail)	4	



- 
 Gentle to Skin
- 
 Helps Protects Against Germ
- 
 Absorbs And Cushions
- 
 Hospital Grade

SOFT & FLEXIBLE Dressing+Pad

For cuts, grazes and covering sutures

This SOFT & FLEXIBLE DRESSING+PAD is really comfortable and conforms to your skin as it is made from a soft medical grade fabric. It has a specially designed absorbent pad which has a cushioning effect to protect the wound. The pad absorbs excess wound fluid which is important for wound healing. The nature of this dressing makes it ideal for most areas of the body, particularly near joints or areas of high movement.

BPOSITIVE ISLAND DRs			
Code	Description / Sizes	PCS/BX	BX/CRT
010604BA	ISLAND DRs 4 X 5CM	50	36
010605BA	ISLAND DRs 5 X 7CM	50	32
010606BA	ISLAND DRs 6 X 8CM	50	24
010608BA	ISLAND DRs 8 X 10CM	25	48
010610BA	ISLAND DRs 10 X 12CM	20	36
010620BA	ISLAND DRs 10 X 20CM	20	18
010605RA	ISLAND DRs 5 X 7CM (Retail)	5	200
010610RA	ISLAND DRs 10 X 12CM (Retail)	5	120
010604RA	Boo Boo Special Dressing 4 x 5cm (PK5) (retail)	5	180



- 
 Gentle to Skin
- 
 Soft and Conformable
- 
 Hospital Grade

SOFT & FLEXIBLE Roll

For retention of dressings and strapping

This SOFT & FLEXIBLE ROLL is really comfortable and conforms to your skin as it is made from a soft medical grade fabric. The widthway stretch of the fabric makes it ideal for most areas of the body, particularly near joints or areas of high movement. Commonly used by Physiotherapist and Podiatrists.

BPOSITIVE FIXATION ROLL			
Code	Description / Sizes	PCS/BX	BX/CRT
010703BA	Soft & Flexible Fixation Roll 2.5CM x 10M	1	128
010705BA	Soft & Flexible Fixation Roll 5CM x 10M	1	48
010705SA	Soft & Flexible Fixation Roll 5CM x 10M (Tan)	1	48
010710BA	Soft & Flexible Fixation Roll 10CM x 10M	1	24
010710SA	Soft & Flexible Fixation Roll 10CM x 10M (Tan)	1	60
010715BA	Soft & Flexible Fixation Roll 15CM x 10M	1	24
010710ZA	Soft & Flexible Fixation Roll 10CM x 2M (Tan)	1	24



Waterproof



Helps Protect Against Germs



Hospital Grade



Waterproof



Helps Protect Against Germs



Hospital Grade



Waterproof



Helps Protect Against Germs



Absorbs And Cushions



Hospital Grade

WATER & BACTERIA PROOF Film Dressing

For cuts, stitches and grazes

This WATER & BACTERIA PROOF Film Dressing is good for minor wounds which are nearly healed and you just wish to give it some extra protection from germs and/or water. Can be used with our Absorbent Dressing to make a customised dressing. Sterile.

BPOSITIVE W/PROOF ISLAND DRS

Code	Description / Sizes	PCS/BX	BX/CRT
020406BA	W/PROOF FILM DRS 6 X 7CM	100	40
020410BA	W/PROOF FILM DRS 10 X 12CM	50	24
020406RA	W/PROOF FILM DRS 6 X 7CM (Retail)	5	120
020406SA	W/PROOF IV U-Shape Film Dressing 6 X 7CM	100	40
020410SA	W/PROOF IV U-Shape Film Dressing 10 X 12CM	50	40
020407SA	W/PROOF IV U-Shape Film Dressing 6.5 X 7CM	100	40

WATER & BACTERIA PROOF Film Roll

For cuts, stitches and grazes

This WATER & BACTERIA PROOF Film Roll helps to protect the wound and reduce the risk of infection. As it is waterproof, you can shower with it.

It is a convenient product to have where you need to cut to a specific length. Non-sterile.

BPOSITIVE W/PROOF FILM ROLL

Code	Description / Sizes	PCS/BX	BX/CRT
020505BA	WATER & BACTERIA-PROOF ROLL 5 X 10M	1	80
020510BA	WATER & BACTERIA-PROOF ROLL 10 X 10M	1	80
020510SA	WATER & BACTERIA-PROOF ROLL 10 x 1M (Retail)	1	300

WATER & BACTERIA PROOF Film Dressing + Pad

For cuts, stitches and grazes

This WATER & BACTERIA PROOF Film Dressing + Pad helps wounds to heal with a specially designed absorbent pad which has a soft cushioning effect to protect the wound. The pad has an extra non-adherent layer which minimise pain when removing the dressing. In addition, our water and bacteria proof protective film helps REDUCE THE RISK OF INFECTION. As it is waterproof, you can shower with it.

BPOSITIVE W/PROOF ISLAND DRS

Code	Description / Sizes	PCS/BX	BX/CRT
020304BA	W/PROOF ISLAND DRS 4 X 5CM	50	40
020305BA	W/PROOF ISLAND DRS 5 X 7CM	50	20
020308BA	W/PROOF ISLAND DRS 8 X 10CM	25	20
020309BA	W/PROOF ISLAND DRS 9 X 20CM	20	20
020312BA	W/PROOF ISLAND DRS 10 X 12CM	25	20
020315BA	W/PROOF ISLAND DRS 10 X 15CM	20	20
020304RA	W/PROOF ISLAND DRS 4 X 5CM (Retail)	5	180
020305RA	W/PROOF ISLAND DRS 5 X 7CM (Retail)	3	180
020308RA	W/PROOF ISLAND DRS 8 X 10CM (Retail)	3	80

RETAIL WOUND CARE PACKS

Our Absorbent Non-Adherent Dressing, Soft & Flexible Dressing and Water & Bacteria Proof Dressing comes in a convenient retail pack for your patients.



*All prices listed above are before GST

Code	Description / Sizes	PCS/BX	BX/CRT
BPOSITIVE N/ADH ABSORBENT			
010105RA	N/ADH ABSORBENT 5 X 5CM (Retail)	4	270
010110RA	N/ADH ABSORBENT 10 X 10 (Retail)	3	120
BPOSITIVE ISLAND DRS			
010605RA	ISLAND DRS 5 X 7CM (Retail)	5	200
010110RA	ISLAND DRS 10 X 12CM (Retail)	5	120
010604RA	Boo Boo Special Dressing 4 x 5CM (PK5) (retail)	5	180
BPOSITIVE W/PROOF ISLAND DRS			
020304RA	W/PROOF ISLAND DRS 4 X 5CM (Retail)	5	180
020305RA	W/PROOF ISLAND DRS 5 X 7CM (Retail)	3	180
020308RA	W/PROOF ISLAND DRS 8 X 10CM (Retail)	3	80
BPOSITIVE W/PROOF DRS			
020306RA	W/PROOF ISLAND DRS 6 X 7CM (Retail)	5	120
BPOSITIVE W/PROOF FILM ROLL			
020510SA	WATER & BACTERIA-PROOF ROLL 10 x 1M (Retail)	1	120
BPOSITIVE SAH HYDROCOLLOID DRS			
061905RA	SAH HYDROCOLLOID D 5CM x 5CM (retail)	5	
061908RA	SAH HYDROCOLLOID D 7.5CM x 7.5CM (retail)	5	
061910RA	SAH HYDROCOLLOID D 10CM x 10CM (retail)	3	
061904RA	SAH HYDROCOLLOID OVAL D 4.4 x 6.9CM (retail)	5	
BPOSITIVE SW CLOSURE STRIPS			
022206RA	SW CLOSURE STRIPS 6MM x 38MM (Retail)	4	

We are here to make a positive impact together, through our products and by giving to worthy causes via **buy1help1**. To achieve this, we have been careful to ensure that our products are of high quality and offer value for money.



BPositive Pte Ltd
Singapore Business Federation Center
160 Robinson Road, #14-04
068914



www.bpositivehealthcare.com